

Daily Chamber Checklist (Perry-Computer)

Please use a separate checklist for each chamber. Initial each line as completed.

| | MON | TUE | WED | THU | FRI |
|---|-----|-----|-----|-----|-----|
| Week Of: | | | | | |
| Initials of person completing checklist | | | | | |
| Chamber # _____ | | | | | |
| Cycle # | | | | | |
| Oxygen zone valves outside & inside the chamber room is on | | | | | |
| Oxygen supply pressure reading @ alarm panel & chamber console: (50-90 psi) | | | | | |
| Air supply pressure reading @ cylinder gauge: (60 psi) if using demand valve, <u>if using free flow system</u> cylinder gauge set at 25 psi | | | | | |
| Chamber covers completely removed and stored in a professional manner | | | | | |
| Chamber ground is connected. Multi meter reading is less than 1 ohm | | | | | |
| Chamber supply/vent hoses without kinks or damage | | | | | |
| Inspect chamber hull for scratches and/or crazing | | | | | |
| Outer ring of the manual override is in the out position | | | | | |
| Electrical Emergency Stop Switch is in the outward position | | | | | |
| Turn the COMMS panel power switch to on and check that the system operates properly | | | | | |
| Flip System START/STOP switch to START . Confirm that indicator eye is showing green. Turn switch off | | | | | |
| Press the (console) main power switch to the depressed position and verify that the red "Power On" indicator light is illuminated | | | | | |
| Press the front panel "Reset" switch and verify the green reset indicator is illuminated | | | | | |
| System "ON" Indicator – Glows green when system is switched to the ON position | | | | | |
| Log on to the chamber computer system | | | | | |
| With door open, rotate swing arm lever to closed position and turn master valve on/check air flow and function of intercom | | | | | |
| Inspect chamber door gasket for damage | | | | | |
| Inspect chamber controls for damage or loose knobs | | | | | |
| Inspect chamber interior and exterior for cleanliness | | | | | |
| Air break equipment disinfected | | | | | |
| Air break cylinder pressure level checked (change below 500psi) | | | | | |

| | | | | | |
|---|--|--|--|--|--|
| Get a treatment started: 1. Select the patient 2. Select the patient's treatment profile 3. Press the RUN TREATMENT button. This will take you to the next screen RUNNING A TREATMENT | | | | | |
| Observation of the control panel (LC) display visually confirms that the treatment profile is actually transpiring | | | | | |
| Adjust ventilation rates (on Chamber Console) from 125 to 385 liters per minute (lpm). Minimum chamber vent setting is 125 LPM | | | | | |
| SHUTDOWN CHECKLIST | | | | | |
| Turn Oxygen and Air supply to chamber off | | | | | |
| Flip System ON/OFF switch to OFF. Confirm that indicator eye is showing NOT showing green | | | | | |
| Turn the COMMS panel power switch off | | | | | |
| Log off the computer | | | | | |
| Press the power button in for the off position and make sure that the power button light goes out | | | | | |
| Cover chamber to protect the acrylic pressure cylinder | | | | | |

Person Completing: _____